


NHS FACT SHEET FOR NHS PROFESSIONALS

The PelvicToner™ progressive resistance vaginal exerciser

The PelvicToner™ offers much improved patient outcomes and significant cost benefits to the NHS in all those situations where pelvic floor exercises are recommended, e.g. stress incontinence, prolapse and sexual dysfunction following childbirth and menopause.

It offers the quickest, simplest and cheapest way to achieve and maintain a stronger pelvic floor and better bladder control, because it encourages subjects to undertake and maintain a pelvic floor exercise programme in the most effective way possible.

The pro-active use of an effective exercise regimen should resolve over 80% of cases of stress incontinence in 8 weeks and obviate the need for pharmaceutical and surgical intervention in the majority of cases. Continued, less frequent use should prevent recurrence.



Applications

The PelvicToner is designed to maximise the effectiveness of pelvic floor muscle exercise and as a training device to ensure a correct exercise regime and maximum compliance.

The PelvicToner is recommended in all the most common uses of pelvic floor exercise, specifically:

- *post-natal* pelvic floor rehabilitation
- primary treatment of *urinary stress incontinence* in women
- treatment of mild-to-moderate cases of *vaginal prolapse* or to supplement other treatments for prolapse that are more serious
- *non-surgical tightening* of the vaginal muscles in cases of sexual dysfunction

Clinical Justification

In clinical trials¹ the PelvicToner was assessed as 'not inferior' to the NICE recommended treatment of Supervised Pelvic Floor Muscle Training (PFMT), but significant additional benefits were noted. Specifically, the PelvicToner:

- aids the subject to confidently identify the correct muscles - the pubo-coccygeal muscle
- ensures an effective exercise regime by providing a resistance that can be varied according to ability
- provides positive feedback during exercise activity and reinforces progressive improvements
- increases compliance through rapid and measurable improvement

In addition, at £15.00 it is less than 1/10th of the cost of the NICE recommended treatment².

Available on prescription

The PelvicToner is available on NHS Prescription and is listed within the Drug Tariff Part IXA-Appliances in the category of Pelvic Toning Devices. NHS Price £15.00

How to find the PelvicToner on EMISS

The PelvicToner can be easily prescribed on EMIS Web, just enter "pelvic" or "pelvic toner" on the medication screen and select PelvicToner from the results.

Distributors: AAH, Alliance Healthcare / Unichem, LexonUK, Phoenix Healthcare Distribution
PIP Code 359-4520

More information

Visit www.pelvictoner.info or contact PT Direct Ltd

¹ Bristol Urological Institute published in *British Journal of Urology International* January 2010)

² NICE CG40 October 2006 pp 161-162 PFMT + Biofeedback costs £168



The PelvicToner is registered as a Class 1 medical device (vaginal trainer) under the provisions of the Medical Devices Directive 93/42/EEC

